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1. Carraro P *et al.* Errors in a stat laboratory. Types and frequencies 10 years later. Clin Chem 2007; 53,7:1338-42.
2. Wennecke G. Useful tips to avoid preanalytical errors in blood gas testing: electrolytes. www.acutecaretesting.org Oct 2003.

Copyright © 2024 Radiometer Medical ApS, Denmark.

Contents may be freely reproduced if the source is acknowledged.

Printed in Denmark by Radiometer Medical ApS, 2700 Brønshøj,
2024.

IVD In vitro diagnostic medical device.

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


1. World Health Organization Regional Office for Europe (2010):WHO guidelines on drawing blood. Best practices in phlebotomy.Geneva, Switzerland: Safe Injection Global Network, World Health Organization, 2010.
2. CLSI.Blood gas and pH analysis and related measurements; approved guideline – second edition.CLSI document C46-A2.2009.

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
1. CLSI. Blood gas and pH analysis and related measurements; approved guideline – second edition. CLSI document C46-A2.2009.
2. World Health Organization Regional Office for Europe (2010): WHO guidelines on drawing blood. Best practices in phlebotomy. Geneva, Switzerland: Safe Injection Global Network, World Health Organization, 2010.
3. Dukic L *et al.* Blood gas testing and related measurements: National recommendations on behalf of the Croatian Society of Medical Biochemistry and Laboratory Medicine. *Biochemia Medica* 2016; 26, 3:318–36.
4. Grenache DG, Parker C. Integrated and automatic mixing of whole blood: An evaluation of a novel blood gas analyzer. *Clin Chim Acta* [Internet]. 2007 Jan 1 [cited 2019 Feb 13];375(1–2):153–7.
5. IFCC guideline Nov 2018 - A primer of Point of Care Blood Gas Testing for laboratorians.

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1. World Health Organization Regional Office for Europe (2010): WHO guidelines on drawing blood. Best practices in phlebotomy. Geneva, Switzerland: Safe Injection Global Network, World Health Organization. CLSI. Blood gas and pH analysis and related measurements; approved guideline – second edition. CLSI document C46-A2. 2009.

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
1. Kahn S. Specimen mislabeling: A significant and costly cause of potentially serious medical errors. www.acutecaretesting.org Apr 2005.
2. Joint Commission on Accreditation of Healthcare Organization. National Patient Goals. www.jointcommission.org Accessed Jun 2013.

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
1. Higgins, C. (2016): Useful tips to avoid preanalytical errors in blood gas testing: pH, pCO_2 and pO_2 . In www.acutecaretesting.org
2. Clinical and Laboratory Standards Institute (CLSI). Blood gas and pH Analysis and Related Measurements; Approved Guideline - second edition. CLSI document C46-A2 (ISBN 1-56238-694-8).
3. WHO guidelines on drawing blood: best practices in phlebotomy. World Health Organization 2010. ISBN 978 92 4 159922 1.

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- Clinical and Laboratory Standards Institute (CLSI). Blood gas and pH Analysis and Related Measurements; Approved Guideline - second edition. CLSI document C46-A2 (ISBN 1-56238-694-8).
- Higgins C. The use of heparin in preparing samples for blood gas analysis. www.acutecaretesting.org Apr 2007.
- Wennecke G. Useful tips to avoid preanalytical errors in blood gas testing: electrolytes. www.bloodgas.org 2003

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
Dukic L *et al.* Blood gas testing and related measurements: National recommendations on behalf of the Croatian Society of Medical Biochemistry and Laboratory Medicine. *Biochemia Medica* 2016; 26, 3:318–36.

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1. Higgins C. The use of heparin in preparing samples for blood gas analysis. www.acutecaretesting.org Apr 2007.
2. Clinical and Laboratory Standards Institute (CLSI). Blood gas and pH Analysis and Related Measurements; Approved Guideline - second edition. CLSI document C46-A2 (ISBN 1-56238-694-8).

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- Toffaletti J. Effect of small air bubbles on changes in blood pO_2 and blood gas parameters: calculated vs. measured effects. www.acutecaretesting.org Jul 2012.
- Higgins C. Pneumatic tube transport of samples for blood gas analysis. www.acutecaretesting.org Jan 2005.

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

1. CLSI. Blood gas pH analysis and related measurements; approved guideline – second edition. CLSI document C46-A2. Clinical and Laboratory Standards Institute, 950 West Valley Road, Suite 2500, Wayne, Pennsylvania 19087 USA, 2009.
2. Narayanan S. Preanalytical issues related to blood sample mixing. www.acutecaretesting.org Oct 2005.
3. Benoit M *et al.* Evaluation and advantages of an automatic magnetic mixing of syringes integrated to a whole blood gas analyser. *Scand J Clin Lab Invest* 2009; 69,5: 628-32.
4. Grenache D *et al.* Integrated and automatic mixing of whole blood: an evaluation of a novel blood gas analyzer. *Clin Chem Acta* 2007; 375:153-57.

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
- Narayanan S. Preanalytical issues related to blood sample mixing. www.acutecaretesting.org Oct 2005.
- Grenache D *et al.* Integrated and automatic mixing of whole blood: an evaluation of a novel blood gas analyzer. *Clin Chem Acta* 2007; 375:153-57.

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1. Wennecke G. Useful tips to avoid preanalytical errors in blood gas testing: electrolytes. www.acutecaretesting.org Oct 2003.
2. Clinical and Laboratory Standards Institute (CLSI). Blood gas and pH Analysis and Related Measurements; Approved Guideline - second edition. CLSI document C46-A2 (ISBN 1-56238-694-8).

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1. Skurup A. Storage recommendations for blood gas samples. Bulletin no. 31. Radiometer Medical ApS publication 2005.
2. Clinical and Laboratory Standards Institute (CLSI). Blood gas and pH Analysis and Related Measurements; Approved Guideline - second edition. CLSI document C46-A2 (ISBN 1-56238-694-8).

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

1. World Health Organization Regional Office for Europe (2010): WHO guidelines on drawing blood. Best practices in phlebotomy. Geneva, Switzerland: Safe Injection Global Network, World Health Organization, 2010.
2. CLSI. Blood gas and pH analysis and related measurements; approved guideline – second edition. CLSI document C46-A2.2009.

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1. CLSI. Blood gas and pH analysis and related measurements; approved guideline – second edition. CLSI document C46-A2.2009.
2. World Health Organization Regional Office for Europe (2010): WHO guidelines on drawing blood. Best practices in phlebotomy. Geneva, Switzerland: Safe Injection Global Network, World Health Organization, 2010.
3. Dukic L *et al*. Blood gas testing and related measurements: National recommendations on behalf of the Croatian Society of Medical Biochemistry and Laboratory Medicine. *Biochemia Medica* 2016; 26, 3:318–36.
4. Grenache DG, Parker C. Integrated and automatic mixing of whole blood: An evaluation of a novel blood gas analyzer. *Clin Chim Acta* [Internet]. 2007 Jan 1 [cited 2019 Feb 13]; 375(1–2): 153–7.

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
Dukic *et al.* Blood gas testing and related measurements: National recommendations on behalf of the Croatian Society of Medical Biochemistry and Laboratory Medicine. *Biochemia Medica* 2016; 26, 3:318–36.

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1. Higgins C. The use of heparin in preparing samples for blood gas analysis. www.acutecaretesting.org Apr 2007.
2. Clinical and Laboratory Standards Institute (CLSI). Blood gas and pH Analysis and Related Measurements; Approved Guideline - second edition. CLSI document C46-A2 (ISBN 1-56238-694-8).
3. half of the Croatian Society of Medical Biochemistry and Laboratory Medicine. *Biochemia Medica* 2016; 26, 3:318–36.

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

1. CLSI. Blood gas pH analysis and related measurements; approved guideline – second edition. CLSI document C46-A2. Clinical and Laboratory Standards Institute, 950 West Valley Road, Suite 2500, Wayne, Pennsylvania 19087 USA, 2009.
2. Narayanan S. Preanalytical issues related to blood sample mixing. www.acutecaretesting.org Oct 2005.
3. Benoit M *et al.* Evaluation and advantages of an automatic magnetic mixing of syringes integrated to a whole blood gas analyser. *Scand J Clin Lab Invest* 2009; 69,5:628-32.
4. Grenache D *et al.* Integrated and automatic mixing of whole blood: an evaluation of a novel blood gas analyzer. *Clin Chem Acta* 2007; 375:153-57.

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
Narayanan S. Preanalytical issues related to blood sample mixing.
www.acutecaretesting.org Oct 2005.

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1. Wennecke G. Useful tips to avoid preanalytical errors in blood gas testing: electrolytes. www.acutecaretesting.org Oct 2003.
2. Clinical and Laboratory Standards Institute (CLSI). Blood gas and pH Analysis and Related Measurements; Approved Guideline - second edition. CLSI document C46-A2 (ISBN 1-56238-6

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1. Dukic L. *et al*, Blood gas testing and related measurements: National recommendations on behalf of the croatian society of medical biochemistry and laboratory medicine, Juli 2016
2. World Health Organization Regional Office for Europe (2010): WHO guidelines on drawing blood. Best practices in phlebotomy. Geneva, Switzerland: Safe Injection Global Network, World Health Organization, 2010.

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1. Kahn S. Specimen mislabeling: A significant and costly cause of potentially serious medical errors. www.acutecaretesting.org Apr 2005.
2. The Joint Commission. National Patient Safety Goals Effective January 2018

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1. WHO guidelines on drawing blood: Best Practices in Phlebotomy; WHO 2010; ISBN 978 92 4 159922 1
2. Clinical and Laboratory Standards Institute (CLSI). Blood gas and pH Analysis and Related Measurements; Approved Guideline - second edition. CLSI document C46-A2 (ISBN 1-56238-694-8).
3. Dukic, L.; Kopčinović, L. M.; Dorotić, A.; Baršić, I. (2016): Blood gas testing and related measurements: National recommendations on behalf of the Croatian Society of Medical Biochemistry and Laboratory Medicine. In *Biochem Med* 26 (3), pp. 318–336. DOI: 10.11613/BM.2016.036.
4. Higgins C. Capillary blood gases – to arterialize or not. www.acutecaretesting.org Jul 2008.

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Higgins C. The use of heparin in preparing samples for blood gas analysis.
www.acutecaretesting.org Apr 2007.

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1. Clinical and Laboratory Standards Institute (CLSI). Blood gas and pH Analysis and Related Measurements; Approved Guideline - second edition. CLSI document C46-A2 (ISBN 1-56238-694-8).
2. Dukic, L.; Kopčinović, L. M.; Dorotić, A.; Baršić, I. (2016): Blood gas testing and related measurements: National recommendations on behalf of the Croatian Society of Medical Biochemistry and Laboratory Medicine. In Biochem Med 26 (3), pp. 318–336. DOI: 10.11613/BM.2016.036.
3. Higgins C. The use of heparin in preparing samples for blood gas analysis. www.acutecaretesting.org Apr 2007.
4. Wennecke G. Useful tips to avoid preanalytical errors in blood gas testing: electrolytes. www.bloodgas.org 2003

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1. Higgins C. The use of heparin in preparing samples for blood gas analysis. www.acutecaretesting.org Apr 2007.
2. Clinical and Laboratory Standards Institute (CLSI). Blood gas and pH Analysis and Related Measurements; Approved Guideline - second edition. CLSI document C46-A2 (ISBN 1-56238-694-8).
3. Dukic L *et al.* Blood gas testing and related measurements: National recommendations on behalf of the Croatian Society of Medical Biochemistry and Laboratory Medicine. *Biochemia Medica* 2016; 26, 3:318-36.

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
1. Dukic L *et al.* Blood gas testing and related measurements: National recommendations on behalf of the Croatian Society of Medical Biochemistry and Laboratory Medicine. *Biochem Med (Zagreb)* 2016; 26,3:318-36.
2. Toffaletti J. Effect of small air bubbles on changes in blood pO_2 and blood gas parameters: calculated vs. measured effects. www.acutecaretesting.org Jul 2012.
3. Higgins C. Pneumatic tube transport of samples for blood gas analysis. www.acutecaretesting.org Jan 2005.
4. Clinical and Laboratory Standards Institute (CLSI). Blood gas and pH Analysis and Related Measurements; Approved Guideline - second edition. CLSI document C46-A2 (ISBN 1-56238-694-8).

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
1. Narayanan, S. Preanalytical issues related to blood sample mixing. www.acutecaretesting.org Oct 2005.
2. Clinical and Laboratory Standards Institute (CLSI). Blood gas and pH Analysis and Related Measurements; Approved Guideline - second edition. CLSI document C46-A2 (ISBN 1-56238-694-8).

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
1. Wennecke G. Useful tips to avoid preanalytical errors in blood gas testing: electrolytes. www.acutecaretesting.org Oct 2003.
2. Clinical and Laboratory Standards Institute (CLSI). Blood gas and pH Analysis and Related Measurements; Approved Guideline - second edition. CLSI document C46-A2 (ISBN 1-56238-694-8).
3. Dukic L *et al.* Blood gas testing and related measurements: National recommendations on behalf of the Croatian Society of Medical Biochemistry and Laboratory Medicine. *Biochem Med (Zagreb)* 2016; 26,3:318-36. この記事は acutecaretesting.org に転載記事としても掲載されています（2017年4月公開）。

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Radiometer Bulletin No. 31. Storage recommendations for blood gas samples
by Anne Skurup. 2005

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1. Clinical and Laboratory Standards Institute (CLSI). Blood gas and pH Analysis and Related Measurements; Approved Guideline - second edition. CLSI document C46-A2 (ISBN 1-56238-694-8).
2. Dukic L *et al.* Blood gas testing and related measurements: National recommendations on behalf of the Croatian Society of Medical Biochemistry and Laboratory Medicine. *Biochem Med (Zagreb)* 2016; 26,3:318–36. この記事は acute-care-testing.org に転載記事としても掲載されています (2017年4月公開)。

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1. Radiometer Bulletin No. 31.Storage recommendations for blood gas samples by Anne Skurup.2005
2. Clinical and Laboratory Standards Institute (CLSI).Blood gas and pH Analysis and Related Measurements; Approved Guideline - second edition. CLSI document C46-A2 (ISBN 1-56238-694-8).
3. Dukic L *et al.* Blood gas testing and related measurements: National recommendations on behalf of the Croatian Society of Medical Biochemistry and Laboratory Medicine. *Biochem Med (Zagreb)* 2016; 26,3:318–36. この記事は acutecaretesting.org に転載記事としても掲載されています（2017年4月公開）。

[illegible][illegible]

[illegible][illegible]

Whatever comes next, we make sure life comes first.

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