**SEPSIS** DETECTION

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
図10: 以下から引用。Kim SJ, et al. Procalcitonin as a diagnostic marker for sepsis/septic shock in the emergency department; a study based on Sepsis-3 definition. American Journal of Emergency Medicine 2018.


図11: 以下から引用。Kim SJ, et al. Procalcitonin as a diagnostic marker for sepsis/septic shock in the emergency department; a study based on Sepsis-3 definition. American Journal of Emergency Medicine 2018.

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1. Singer M et al. The Third International Consensus Definitions for Sepsis and Septic Shock (Sepsis-3). *JAMA* 2016; 315: 801-10
2. Rhodes A et al. Surviving Sepsis Campaign: International Guidelines for Management of Sepsis and Septic Shock: 2016. *Intensive Care Med* 2017; 43.
3. Hotchkiss RS et al. Sepsis and septic shock. *Nat Rev Dis Primers* 2016; 3 0: 2.
4. Hall MJ, Williams SJ, DeFrances CJ, Golosinsky A. Inpatient care for septicemia or sepsis: A challenge for patients and hospitals. Available at: <http://www.cdc.gov/nchs/data/databriefs/db62.pdf> Last assessed June 2020.
5. Cohen J, Vincent JL, Adhikari NK et al. Sepsis: a roadmap for future research. *Lancet Infect Dis* 2015; 15.
6. Vincent JL, De Backer D. Circulatory shock. *N Engl J Med* 2013; 369:
7. Wang HE, Shapiro NI, Angus DC, Yealy DM. National estimates of severe sepsis in United States emergency departments. *Crit Care Med* 2007; 35:
8. Kumar A, Roberts D, Wood KE et al. Duration of hypotension before initiation of effective antimicrobial therapy is the critical determinant of survival in human septic shock. *Crit Care Med* 2006; 34:
9. Sterling SA, Ryan Miller W, Pryor J et al. The impact of timing of antibiotics on outcomes in severe sepsis and septic shock: A systematic review and metaanalysis. *Crit Care Med* 2015; 43:
10. Ferrer R, Martin-Loeches I, Phillips G et al. Empiric antibiotic treatment reduces mortality in severe sepsis and septic shock from the first hour: results from a guideline-based performance improvement program. *Crit Care Med* 2014; 42:
11. Gaieski DF, Mikkelsen ME, Band RA et al. Impact of time to antibiotics on survival in patients with severe sepsis or septic shock in whom early goal-directed therapy was initiated in the emergency department. *Crit Care Med* 2010; 38:
12. Lee SH, Chan RC, Wu JY et al. Diagnostic value of procalcitonin for bacterial infection in elderly patients - a systemic review and meta-analysis. *Int J Clin Pract* 2013; 67.
13. Meisner M. Update on procalcitonin measurements. *Ann Lab Med* 2014; 34.
14. Meisner, M et al. Procalcitonin--influence of temperature, storage, anticoagulation and arterial or venous asservation of blood samples on procalcitonin concentrations. *Eur J Clin Chem Clin Biochem* 1997; 35, 8.
15. Sakr Y, Sponholz C, Tuche F, Brunkhorst F, Reinhart K. The role of procalcitonin in febrile neutropenic patients: review of the literature. *Infection* 2008; 36.
16. Meisner M. Procalcitonin: Erfahrungen mit einer neuen Meßgröße für bakterielle Infektionen und systemische Inflammation. *J Lab Med* 1999, 23 (5); 263-272
17. Tang H, Huang T, Jing J, Shen H, Cui W. Effect of procalcitonin-guided treatment in patients with infections: a systematic review and meta-analysis. *Infection* 2009; 37.
18. Jones AE, Fiechtl JF, Brown MD, Ballew JJ, Kline JA. Procalcitonin test in the diagnosis of bacteremia: a meta-analysis. *Ann Emerg Med* 2007; 50.
19. Becker KL, Snider R, Nylen ES. Procalcitonin in sepsis and systemic inflammation: a harmful biomarker and a therapeutic target. *Br J Pharmacol* 2010; 159.
20. Freund Y et al. Serum lactate and procalcitonin measurements in emergency room for the diagnosis and risk-stratification of patients with suspected infection. *Biomarkers* 2012; 17.
21. Riedel S, Melendez JH, An AT, Rosenbaum JE, Zenilman JM. Procalcitonin as a marker for the detection of bacteremia and sepsis in the emergency department. *Am J Clin Pathol* 2011; 135:
22. Balk RA et al. Effect of procalcitonin testing on health-care utilization and costs in critically ill patients in the United States. *Chest* 2017; 151:
23. Schuetz P. et al. Procalcitonin (PCT)-guided antibiotic stewardship: an international experts consensus on optimized clinical use. *Clin Chem Lab Med* 2019; 57, 9.
24. Daniels, R.; Nutbeam, T.; McNamara, G.; Galvin, C. The sepsis six and the severe sepsis resuscitation bundle: A prospective observational cohort study. *Emerg. Med. J.* 2011, 28.
25. National Institute for Health and Care Excellence. Sepsis: Recognition, Diagnosis and Early Management (NICE Guideline NG51). 2016. Available online: <https://www.nice.org.uk/guidance/ng51> (last accessed June 2020)
26. Levy MM, Evans LE, Rhodes A. The Surviving Sepsis Campaign Bundle: 2018 update. *Intensive Care Med* (2018) 44:925-928
27. Opal S, Wittebole X. Biomarkers of Infection and Sepsis. *Crit Care Clin* 2020; 36.
28. Broyles M. Impact of Procalcitonin-Guided Antibiotic Management on Antibiotic Exposure and Outcomes: Real-world Evidence. *Open Forum Infectious Diseases* 2017
29. Teggert A, Datta H, Ali Z. Biomarkers for Point-of-Care Diagnosis of Sepsis. *Micromachines* 2020, 11, 286.
30. Kondo, Y.; Umemura, Y.; Hayashida, K.; Hara, Y.; Aihara, M.; Yamakawa, K. Diagnostic value of procalcitonin and presepsin for sepsis in critically ill adult patients: A systematic review and metaanalysis. *J. Intensive Care* 2019; 7.
31. Christ-Crain M, Stolz D, Bingisser R et al. Procalcitonin guidance of antibiotic therapy in community-acquired pneumonia: a randomized trial. *AM J Respir Crit Care Med* 2006; 174:
32. Schuetz P, Briel M, Christ-Crain M et al. Procalcitonin to guide initiation and duration of antibiotic treatment in acute respiratory infections: an individual patient data meta-analysis. *Clin Infect Dis* 2012; 55.
33. Christ-Crain M, Jaccard-Stolz D, Bingisser R et al. Effect of procalcitonin-guided treatment on antibiotic use and outcome in lower respiratory tract infections: cluster-randomised, singleblinded intervention trial. *Lancet* 2004; 363.
34. Bouadma L, Luyt C-E, Tubach F et al. Use of procalcitonin to reduce patients' exposure to antibiotics in intensive care units (PRORATA trial): a multicentre randomised controlled trial. *Lancet* 2010; 375.
35. Gould IM. Controversies in infection: infection control or antibiotic stewardship to control healthcare-acquired infection? *J Hosp Infect* 2009; 73.
36. Maisel A, Neath S-X, Landsberg J et al. Use of procalcitonin for the diagnosis of pneumonia in patients presenting with a chief complaint of dyspnoea: results from the BACH (Biomarkers in Acute Heart Failure) trial. *Eur J Heart Failure* 2012; 14.
37. Kim SJ, et al. Procalcitonin as a diagnostic marker for sepsis/septic shock in the emergency department; a study based on Sepsis-3 definition. *Am J Emerg Med* 2018

## References continued:

38. Ferrer R, Artigas A, Suarez D, Palencia E, Levy MM, Arenzana A, et al. Effectiveness of treatments for severe sepsis: a prospective, multicenter, observational study. *Am J Respir Crit Care Med* 2009;180:861–6.

39. Cervellin G, Schuetz P, Lippi G. Toward a holistic approach for diagnosing sepsis in the emergency department. Chapter 5 in *Advances in Clinical Chemistry* 2019; 92

40. Haydar S, Spanier M, Weems P, Wood S, Strout T. Comparison of QSOFA score and SIRS criteria as screening mechanisms for emergency department sepsis. *Am J Emerg Med* 2017;35:1730–3.

41. Liu SW, Thomas SH, Gordon JA, Hamedani AG, Weissman JS. A pilot study examining undesirable events among emergency department-boarded patients awaiting inpatient beds. *Ann Emerg Med*

## Whatever comes next, we make sure life comes first

Radiometer products and solutions are used in hospitals, clinics, and laboratories in over 130 countries, providing information on critical parameters in acute care diagnostics. Through connected solutions, expert knowledge, and trusted partnership, we help health care professionals make diagnostic decisions to improve patient care.

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© Radiometer Medical ApS, 2700 Brønshøj, Denmark, 2020. All Rights Reserved. MAPSSS 000412 R1.